# Selectra – Accessory Kit

---

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### General considerations

#### Target group

This manual is intended for cardiologists, electrophysiologists and surgeons who are familiar with and have practical experience with cardiac resynchronization therapy and implantation of leads in the coronary venous system.

#### Description

BIOTRONIK's Selectra CS lead introducer system is a combination of guiding catheters and implantation accessories used to facilitate access to the coronary venous system for suitable leads and catheters.

#### Packaging units

It primarily consists of the accessory kit and various, individually available guiding catheters with a nominal inner diameter of 7 F (outer guiding catheters) or 5 F (inner guiding catheters).

#### Advantages for the user

Using the Selectra system has the following advantageous options:

- Quick probing of the coronary sinus
- Easy access to the coronary sinus
- Introduction of contrast media for angiography
- Easy and fast positioning and repositioning of the leads

### Intended use

The Selectra accessory kit is used in conjunction with the CS lead introducer system to facilitate lead implantation in the left side of the heart via the coronary sinus.

### Contraindications

Use of the Selectra CS lead introducer system is contraindicated for the following:

- Patients with an existing or possible occlusion of the coronary vessels or unsuitable anatomy of the coronary veins
- Patients with active systemic infection

### Possible negative side effects

Possible complications are as follows:

- Allergic reactions to contrast media
- Hematoma
- Bleeding
- Infection
- Embolism
- Pneumothorax
- Cardiac tamponade
- Myocardial damage
- Venous or cardiac perforation
- Local tissue response, fibrotic tissue formation
- Damage to the heart valves
- Vascular occlusion
- Chronic nerve damage

### Packaging, sterility, storage, and disposal

#### Sterile container

The Selectra accessory kit is delivered in a sterile container consisting of a blister with sealing paper. This blister is sterily sealed and placed in a sterile bag, so that the outside of the blister is also sterile. Ethylene oxide is used for sterilization purposes.

#### Outer box and labeling

The sterile bag with the blister is packaged in a box that has a quality control seal and a product information label. The label lists the model type, technical data, use by date, and information about sterility and storage of the package and content.

**1** To ensure sterility, the container should be checked for damage prior to opening.

**2** If it is suspected that the sterile packaging has been opened or damaged, return the product to BIOTRONIK.

#### Storage

The permissible storage temperature ranges between 5 and 55°C.

If this temperature range or the use by date are exceeded, the documented properties of the implantation accessories can no longer be warranted. This could result in technical failures.

#### Disposal

Used implantation accessories must be disposed of as medical waste in an environmentally friendly and proper manner.

They do not contain material that requires further provisions.

### Scope of delivery (package contents)


Fig. 1: Accessory kit contents

### Specifications

Amount, designation	Specification
1 slitter tool	For leads with a diameter of between 4.6 and 5.8 F Can cut up to 3 Selectra catheters
1 guide wire (Seldinger wire)	Length: 150 cm Material: stainless steel, teflon coated Diameter: 0.89 mm
4 transvalvular insertion tools (TVI)	Inner diameter: 0.42 mm [1.26 F]
1 syringe	With Luer lock connection, volume: 12 cm <sup>3</sup>
1 torque tool for OTW guide wires	For guide wires with 0.36 mm diameter (for OTW leads)
2 one-way valves	For the side port, bidirectional, with Luer lock internal and external thread. Opened by screwing in syringe
2 one-way stopcocks	For the side port, with Luer lock internal and external thread
2 sealing caps	For the side port, with Luer lock inner thread
1 technical manual	Outside the sterile container

Table 1: Scope of delivery and specifications

## Basic Safety and Handling Instructions

### Not for multiple usage

The Selectra implantation accessory kit is intended for single-use only.

Reuse of implantation accessories can result in infections, embolisms and damage to the product.

### Avoiding thrombosis

To inhibit blood clotting, BIOTRONIK recommends intravenously administering 50-100 IU heparin per kg body weight if the patient is not otherwise anticoagulated. The heparin may be directly administered through the Selectra guiding catheter or a different venous access point. In this case, please rinse with sufficient physiological NaCl solution to guarantee that the heparin completely enters the patient's circulation.

### Damaging the lead

Avoid damaging the lead with sharp-edged instruments - especially with the blade of the slitter tool when removing the catheter.

### Contrast media

Only use water soluble contrast media.

### X-ray monitoring

Generally, perform the described implantation processes only under X-ray monitoring.

### Damaged catheters

Do not use an (inner or outer) catheter showing any signs of damage.

### Rinse catheter thoroughly before use!

The inner and outer Selectra guiding catheters must be rinsed thoroughly before use (using a physiological NaCl solution), otherwise the hydrophilic inner coating of the catheters, which significantly improves their gliding property, is not effective. This will also help to avoid trapping air in the catheters.

### Using the transvalvular insertion tool

To prevent kinking, a transvalvular insertion tool (TVI) has to be used to guide a 0.36 mm OTW lead guide wire through the valve and into the guiding catheter.

## Handling Instructions

**Note:** This technical manual does not claim to explain all the details of the implantation of a CS lead. This technical manual only deals with certain important aspects and examples to illustrate the proper handling of the implantation accessories.

### Opening the packaging

#### In a non-sterile area

Remove the sterile bag from the storage package and open it.

#### In a sterile area

Remove the sterile blister and open it by peeling the sealing paper off in direction of the arrow.

### Preparation

**1** Choose the Selectra guiding catheter or a pairing of inner and outer catheters from the product line that has the best dimensions and geometry for the patient's atrial anatomy.

**2** Please consult the technical manual included with the catheters.

**3** Prior to use, all components must be thoroughly rinsed and, if possible, vented. Repeat the rinsing procedure at regular intervals during the use of the implantation accessories.

**Note:** Use a physiological NaCl solution for rinsing.

**Note:** The connection of a hemostatic valve is not required as all Selectra guiding catheters have a hemostatic valve embedded in the handle.

**4** The side port of the guiding catheter must be closed by fixing either the sealing cap, the one-way stopcock, or the bidirectional check valve to the Luer lock connection. This applies to inner as well as outer catheters.

### Using the transvalvular insertion tool

To prevent kinking, a transvalvular insertion tool (TVI) has to be used to guide a 0.36 mm OTW lead guide wire through the valve and into the guiding catheter.

**5** Insert the associated dilator into the (outside) guiding catheter. This causes the catheter to straighten.

**Note:** Do not attempt to reshape the guiding catheter! It is preformed (except for curve shapes Straight and Straight L) and resumes its curved shape after the dilator has been removed.

### Puncture

**1** Connect a suitable puncturing cannula (inner diameter at least 1 mm = 18 G) to a syringe body. Puncture the selected vein (e.g., subclavian vein, cephalic vein) at a suitable site.

**Caution!** For patients with chronic lung disease, there is an increased risk of pneumothorax!

**2** To check the venous puncture, aspirate blood into the syringe.

**3** Remove the syringe body. Leave the cannula in position.

## Introduce guiding catheter up to the atrium

1 Using X-ray monitoring, introduce the included guide wire through the cannula into the vein and advance it up to the atrium.

2 Remove the cannula.

**Note:** You can also use a peel-away introducer sheath with a dilator. A lead introducer with an inner diameter of at least 9 F is suitable for this purpose.

- With the dilator fully in place, push the introducer sheath over the guide wire into the vessel.
- Then remove the dilator. The introducer sheath remains in the vessel.

3 Introduce the selected outer guiding catheter and dilator through the guide wire. Advance the guiding catheter up to the atrium.

**Caution!** Carefully advance the guiding catheter to prevent damage to the vessel walls!

4 Remove the dilator and the guide wire. The guiding catheter resumes its reshaped form (except for curve shapes Straight and Straight L), thus facilitating access to the coronary sinus.

## Probing of the coronary sinus, introduce catheter into coronary sinus

The decisive step during the implantation of a lead into the coronary venous system is the introduction of the catheter into the coronary sinus. Different approaches to do this are provided by using the Selectra CS lead introducer system and possibly further accessories.

### Connecting syringe with side port (introduction of contrast media)

The side ports of the inner and outer guiding catheters are identical, as is the process of introducing contrast media for each type.

The syringe includes a Luer lock thread, which can be screwed into the Luer lock thread of the bidirectional valve or one-way stopcock.

If a sealing cap has been covering the catheter's side port, this must be replaced by a bidirectional valve or one-way stopcock.

The one-way stopcock can be opened or closed manually at any given time.

The bidirectional valve can be opened by screwing the Luer lock thread of the syringe. It closes automatically when the syringe is removed.

### Probing of the coronary sinus directly with the outer guiding catheter

1 Push the guiding catheter into the coronary sinus through the ostium (except curve shapes Straight and Straight L).

A one-way stopcock or check valve can be connected to the side port of the guiding catheter and used to funnel in contrast media.

### Probing of the coronary sinus with a straight outer catheter and steerable EP catheter

A straight guiding catheter (Straight or Straight L) can only be formed or guided in the coronary sinus using a steerable EP guiding catheter (e.g., ViaCath NG 4/5/5 mm [not available in the U.S]).

1 Introduce the steerable EP catheter into the straight guiding catheter.

2 Probe the coronary sinus with the steerable EP catheter and push the EP catheter into the coronary sinus.

3 Push the straight guiding catheter into the coronary sinus using the steerable EP catheter. A one-way stopcock or check valve can be connected to the side port of the guiding catheter and used to funnel in contrast media.

### Probing the coronary sinus with inner guiding catheters or diagnostic catheters

Probing of the coronary sinus can be facilitated using the Selectra inner guiding catheter or a suitable diagnostic catheter (e.g., HyperFlow MPA2 or AL1).

Selectra inner and outer catheters make up a telescopic system.

Observe the respective manuals for use of a diagnostic catheter.

A one-way stopcock or check valve can be connected to the side port of the inner and outer Selectra guiding catheter and used to funnel in contrast media if necessary.

**Note:** The main part of the inner Selectra guiding catheter also has a Luer lock, so the syringe with a Luer lock thread can be fixed directly to the main port. Therefore, the inner catheter can be steered single-handedly and a contrast medium can be injected at the same time.

1 Advance the probing catheter (diagnostic catheter or Selectra inner catheter) to the atrium through the implanted [outer] guiding catheter.

2 Probe the coronary sinus with the probing catheter and introduce the probing catheter into the coronary sinus.

3 Advance the guiding catheter over the probing catheter into the coronary sinus or further advance both catheters together.

### Probing of a lateral target vein for the lead using the inner catheter

Carefully insert the tip of the inner catheter into the target vein and then advance it. Insert a suitable guide wire for the designated OTW lead [max. diameter 0.36 mm or 0.014 inch] through the inner guiding catheter into the target vein and, if needed, advance it further.

The outer Selectra catheter can now be advanced further via the inner catheter.

## Angiography

To obtain a better overview of the coronary venous system, a coronary angiography may be carried out.

A contrast medium can be injected directly through a lateral inner catheter. However, the contrast medium spreads out better when the venous blood flow is temporarily interrupted using a balloon catheter.

- Remove inner catheter or EP diagnostic catheter.
- Insert a suitable venogram balloon catheter (e.g., Corodyn P1 [not available in the U.S]) into the coronary sinus using the guiding catheter.

Please consult the user manual provided with the venogram balloon catheter.

Keep the interruption of the blood flow as short as possible.

## Transvalvular insertion tool (TVI)

The transvalvular insertion tool (TVI) is guided into the integrated hemostatic valve to allow or facilitate the process of steering 0.36 mm guide wires for OTW leads through the valve.

Use the transvalvular insertion tool for optimal valve effectiveness. Prior to moving the lead forward, remove the transvalvular insertion tool using the peel-away technique.

## Inserting the lead

### Positioning of an OTW lead with a diameter less than 5 F

- An OTW lead with an outer diameter less than 5 F can be steered directly through the inner guiding catheter.

The rest of this process is described in the manual of the respective CS OTW lead.

In order to then remove the inner and outer guiding catheter, the procedure with the slitter tool described below needs to be carried out twice in a row (first for the inner, and then for the outer catheter).

### Positioning of an OTW lead with a diameter greater than 5 F

- Remove the inner guiding catheter, at the same time hold the guide wire in its position.
- Position OTW lead over the guide wire through the outer catheter into the target vein.

The rest of this process is described in the manual of the respective CS OTW lead.

**Note:** The implantation accessories for the Selectra CS lead introducer system can only be used in combination with suitable leads with polyurethane coating or a coating with comparable gliding properties and a diameter of 4.6 to 5.8 F.

## Remove guiding catheter using slitter tool

**Note:** The guiding catheters for BIOTRONIK's Selectra CS lead introducer system must be slit open and removed with a slitter tool.

1 Take the slitter tool in one hand, and the lead with the catheter that is to be removed in the other - see Figure 2.

**Note:** The slitter tool is suitable for both left-handed and right-handed use.


Fig. 2: Holding the slitter tool, clamp not open

Legend:

- 1 Grip plate of the slitter tool
- 2 Blade
- 3 Clamp for the lead

2 To open the clamp, press and hold the release button.

3 With the other hand, insert the lead into the clamp.


Fig. 3: Insert the lead into the open clamp

4 Disengage the release button - the lead is held in the clamp.

5 Position the handle of the guiding catheter on the blade of the slitter tool.

6 Support the hand that is holding the slitter tool - if possible -, holding the catheter handle level with the slitter tool. Pull the guiding catheter against the slitter tool slitting it open.


Fig. 4: Hold the slitter tool level with the catheter handle. Pull the catheter parallel to the edge of the slitter tool against the blade.

If the lead is part of a telescope system consisting of an inner and an outer catheter, first the inner and then the outer catheter is removed using the slitter tool. The slitter tool is suitable for slitting up to three catheters.

**Note:** When using the slitter tool, make sure that the guiding catheter moves straight along the slitter tool (Fig. 4). If the slitting procedure requires an unusual amount of effort or other signs of problems occur while slitting the catheters, resume slitting or reclamp the lead or replace the slitter tool with a new one and then continue slitting.

7 If present, remove the introducer sheath appropriately.

## Product Line for the Selectra CS Lead Introducer System

The Selectra CS lead introducer system consists of the accessory kit described in this technical manual and a selection of compatible guiding catheters with different curve shapes and working lengths. Table 2 shows our product line.

Figure	Item	Order number
<b>Outer guiding catheters for Selectra</b>		
	Selectra Amplatz 6.0-45	375 519
	Selectra Amplatz 6.0-55	375 520
	Selectra BIO 2-45	375 523
	Selectra BIO 2-55	375 524
	Selectra Extended Hook-45	375 527
	Selectra Extended Hook-55	375 528
	Selectra Hook-45	375 529
	Selectra Hook-55	375 530
	Selectra MPEP-45	375 531
	Selectra MPEP-55	375 532
	Selectra Multipurpose Hook-45	375 533
	Selectra Multipurpose Hook-55	375 534
	Selectra Right-45	375 535
	Selectra Right-55	375 536
	Selectra Straight-45	375 537
	Selectra Straight-55	375 521
<b>Inner guiding catheters for Selectra</b>		
	Selectra IC-50-59	392 289
	Selectra IC-50-69	392 290
	Selectra IC-90-59	392 291
	Selectra IC-90-69	392 292
<b>Selectra Accessory Kit</b>		
	Slitter Tool	383 119
<b>Suitable accessories</b>		
Suitable lead introducer set including valve and side port with an internal diameter of at least 9 F		
Steerable EP catheter ViaCath NG 4/5/5mm (not available in the U.S.) (Only for outer catheters)		351 197
Venogram balloon catheter Corodyn P1(not available in the U.S.) (Only for outer catheters)		336 074

Table 2: Available components and accessories for the Selectra CS lead introducer system

### Compatible leads and accessories

The individually available Selectra guiding catheters with a nominal inner diameter of 7 French (outer catheter) and 5 French (inner catheter) can only be used in combination with the accessory kit of the Selectra CS lead introducer system. Outer guiding catheters and implantation accessories of the Selectra CS lead introducer system are suited for leads with polyurethane coatings or coatings with comparable anti-friction properties and a diameter between 4.6 and 5.8 F. Leads with a diameter of less than 5.0 F can be inserted directly through the inner guiding catheter.

Suitable CS guiding catheters and combinations of inner and outer catheters for different lead lengths can be found in the following table.

BIOTRONIK CS leads		Selectra guiding catheter	
Outer diameter	Length	Outer catheter Working length	Inner catheter Working length
≤ 5.8 F	77 cm	45 cm	-
		55 cm	-
<5.0 F	87 cm	45 cm	-
		55 cm	-
	77 cm	45 cm	-
		55 cm	-
	87 cm	45 cm	59 cm
		55 cm	-
	97 cm	45 cm	59 cm
			69 cm
		55 cm	69 cm

Table 3: Compatibility of CS guiding catheters and leads

### Technical Data for the Guiding Catheters

All guiding catheters		Shorter models	Longer models
Material		PEBAX	
X-ray contrast		Barium sulfate in the entire catheter area	
Radiopaque marker at the distal end		Tungsten (black marker zone, 20 mm wide)	
Soft tip length:		1.0 mm (blue)	
Outer guiding catheters			
Working length		45.0 cm	55.0 cm
Total length		48.3 cm	58.3 cm
Inner diameter		2.44 mm [7.3 F]	
Outer diameter		2.91 mm [8.7 F]	
Handle		Ergonomic design, side port with Luer lock connector, main port with integrated hemostatic silicone valve	
Inner guiding catheters			
Working length		59.0 cm	69.0 cm
Total length		62.9 cm	72.9 cm
Inner diameter		1.85 mm [5.55 F]	
Outer diameter		2.31 mm [6.93 F]	
Handle		Ergonomic design, side port with Luer lock connector, main port with integrated hemostatic silicone valve and with Luer lock connector	

<b>Dilator (only for outer catheters)</b>		
Working length	53.5 cm	63.0 cm
Material	PE	
Outer diameter	2.3 mm [6.9 F]	
Inner diameter	1.02 mm [3.06 F]	
Proximal end	Luer lock connection	

Table 4: Technical data of the guiding catheters

### Disclaimer

BIOTRONIK's Selectra CS lead introducer system (including all individual accessories available separately), used in conjunction with devices for implantation in the coronary sinus, has been qualified, manufactured, and tested according to proven and accepted standards and procedures. Nonetheless, the physician must be aware that the accessories may easily be damaged through improper handling or use. Except as set forth in BIOTRONIK's limited warranty, BIOTRONIK makes no expressed or implied warranties for its accessories.

### Legend for the Label

Manufacturing date

Use by

Temperature limit

BIOTRONIK order number

Lot number

Sterilized with ethylene oxide

Do not reutilize

Do not reuse

Non-sterile

Follow the instructions for use

Contents

Do not use if packaging is damaged

CE mark

Inner diameter

Outer diameter

Total length

Working length


Guide wire

Torque tool for OTW guide wires

Dilator

Selectra guiding catheters (various curve shapes)

Selectra transvalvular insertion tool (peel-away introducer TVI)

Slitter tool (Selectra)

Selectra one-way stopcock

Selectra one-way valve

Selectra sealing cap

Syringe

Caution:  
Federal (U.S.A.) law restricts this product to sale by, or on the order of, a physician.

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All rights reserved.
BIOTRONIK, Inc.